




# Viral DNA and RNA Extraction Kit

## User Guide



Version 3.0



In-Vitro Diagnostics / For use with Automatic nucleic acid extractor compatible with  
Viral DNA and RNA Extraction Kit



T108H



**Xi'an Tianlong Science and Technology Co., Ltd.**

No.389, Zhuhong Road, Xi'an, 710018, Shaanxi, P.R. China;

2-3F, No. 2 Building, Yuanzheng Innovation Park, Caotan Ecological Zone, No.1258, Hongye Road,  
Economic & Technological Development Zone, Xi'an, 710018, Shaanxi, P.R. China.



**SUNGO Europe B.V.**

Olympisch Stadion 24, 1076DE Amsterdam, Netherlands.

# Contents

|   |   |
|---|---|
| Intended Use .....                        | 1 |
| Product Performance Indicators .....      | 1 |
| Special Notes .....                       | 1 |
| Testing Principle .....                   | 1 |
| Kit Contents .....                        | 2 |
| Materials Required but not Provided ..... | 2 |
| Warnings and Precautions .....            | 2 |
| Reagent Storage and Handling .....        | 3 |
| Sample Handling and Storage .....         | 3 |
| Operation Guide .....                     | 4 |
| 1. Reagent Preparation .....              | 4 |
| 2. Extraction Steps .....                 | 4 |
| Troubleshooting Guide .....               | 4 |
| Quality Control .....                     | 5 |
| Limitations of Test Methods .....         | 5 |
| Safety Symbols and Signs .....            | 6 |
| Contact Information .....                 | 6 |

## Intended Use

The **Viral DNA and RNA Extraction Kit** is designed to extract viral RNA/DNA from serum, swab etc. The extracted Viral DNA and RNA is of high purity and stability, and can be used in a variety of routine operations, including enzyme digestion, Polymerase Chain Reaction (PCR), DNA library constructions, Southern hybridization and blotting, Next-Generation Sequencing (NGS) and other experiments.

The **Viral DNA and RNA Extraction Kit** is intended to be used by professionals, such as biotechnologists, microbiologists, clinical technicians, and physicians who are trained in molecular and biological techniques.

## Product Performance Indicators

The extraction kit can extract more than 10 IU/mL viral DNA nucleic acid, and more than 30 IU/mL viral RNA nucleic acid. Both the intra and inter-batch variations of kit are less than 5%.

## Special Notes

The extraction kit is particularly used for viral DNA/RNA isolation; therefore, all of experiment supplies, such as pipettes, tubes, tips, must be autoclaved. Operator should wear gloves and masks. Before using the kit, please read the manual and strictly follow the protocol. Clinical samples should be processed on clean bench or in biosafety cabinet.

Do not mix reagents from different batches, and use the kit within expiry date.

Dispose of all samples and reagent materials used in an experiment, and thoroughly clean and disinfect the experimental workbench.

## Testing Principle

During the process of nucleic acid extraction, using silica gel column principle to conduct adsorption and washing of nucleic acid, which can achieve the transfer of nucleic acid and complete the extraction and purification of nucleic acid.




Figure 1. Schematic Diagram of Spin-Column Method

## Kit Contents

| NO. | Components           | Specification | Quantity |
|-----|----------------------|---------------|----------|
| 1   | Lysis Solution       | 16mL          | 1 bottle |
| 2   | Washing A            | 22.4mL        | 1 bottle |
| 3   | Washing B            | 22.4mL        | 1 bottle |
| 4   | Eluent               | 1.6mL         | 1 bottle |
| 5   | Proteinase K         | 14mg          | 1 bottle |
| 6   | Proteinase K Diluent | 0.7mL         | 1 bottle |
| 7   | Adsorption Column    | /             | 32       |
| 8   | Collection Tube      | /             | 64       |
| 9   | Sponge Inner Support | /             | 1        |
| 10  | Packaging Box        | /             | 1        |
| 11  | Instructions for Use | /             | 1 copy   |

## Materials Required but not Provided

When working in a laboratory, make sure to wear a proper lab coat, powder-free disposable gloves and protective goggles. For more information, please consult the Safety Data Sheet (SDS) available from the product supplier.

- Pipettor: 50µL, 200µL, 1000µL
- Tip: 50µL, 200µL, 1000µL
- 1.5mL nuclease-free centrifuge tube
- Vortex mixer
- High-speed centrifuge
- Water bath
- Sample holder
- 75% ethanol
- EtOH (to be ordered separately)

## Warnings and Precautions

**Please be sure to read the precautions before using the kit.**

The extraction kit is particularly used for viral DNA/RNA isolation. Therefore, all of experiment supplies, such as pipettes, tubes, tips, must be processed by autoclave. Operator should wear gloves and masks.

Please read the manual carefully before using the kit, and strictly follow the manual throughout the operation. The clinical samples should be collected on a clean bench or in a bio-safety cabin.

Do not mix reagents from different batches, and use the kit within expiry date.

After the experiment, dispose of all samples and reagent materials used in an experiment, and clean and disinfect the experimental workbench thoroughly .

The ***Viral DNA and RNA Extraction Kit*** is intended for in vitro diagnosis use.

When using kit, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate Material Safety Data Sheets (MSDSs). These documents are available online in a convenient and compact PDF format at <https://www.ug-msds.com/MSDS1>, where the operator can find, view and print the appropriate MSDSs.

 **Caution: Do not add any bleach or acidic solution directly to the pre-filled reagent.**

The pre-filled reagent contains guanidinium salts, which, when combined with bleach can form highly reactive compounds. If any of buffers are spilled, clean immediately with a suitable laboratory detergent and water. If the spilled liquid contains potentially infectious agents, clean the affected area first with laboratory detergent and water first. Then clean with sodium hypochlorite at a concentration of 1% (v/v). The kit comes with the following warnings and precautions.

● **Product contents**

Guanidine hydrochloride, sodium dodecyl sulfate, trihydroxymethyl aminomethane, isopropanol, absolute ethanol.

● **Toxicological information**

**Skin corrosion/irritation**

May irritate the skin.

**Severe eye damage/eye irritation**

May cause irreversible eye damage.

**Respiratory or dermal sensitivity**

No relevant data is available.

**Germ cell mutagenesis**

Do not conform with the classification criteria based on the existing data.

**Carcinogenicity**

Do not conform with the classification criteria based on the existing data.

**Reproductive toxicity**

Do not conform with the classification criteria based on the existing data.

**Specific target organ toxicity - single exposure**

May cause drowsiness or dizziness.

**Specific target organ toxicity - repeated exposure**

Do not conform with the classification criteria based on the existing data.

**Potential health effects**

Inhalation: Avoid inhalation of concentrated vapour. Inhaling a large amount of vapour may cause respiratory irritation. May cause drowsiness or dizziness.

Skin contact: May cause skin irritation.

Eye contact: Liquid contact may cause eye damage.

Ingestion: For any unexpected route of exposure, it may be harmful if ingested.

● **Ecological Information**

Ecotoxicity: Harmful to aquatic life with long-lasting effects.

Mobility: No relevant data is available.

Bioaccumulation potential: No relevant data is available.

Environmentally adverse effects: No relevant data is available.

Other adverse effects: Do not allow the product to enter drains or water sources.

● **First aid measures**

In case of eye contact: Immediately rinse the upper and lower eyelids with plenty of water.

In case of skin contact: Immediately remove contaminated clothing thoroughly, rinse the skin with soap and plenty of water. If irritation persists, immediately contact the nearest doctor/physician.

In case of inhalation: Keep away from exposure and transfer to a place with fresh air.

In case of ingestion: Do not give anything orally to an unconscious person. Rinse mouth thoroughly with water and seek immediate medical attention for symptomatic treatment.



## Reagent Storage and Handling

The **Viral DNA and RNA Extraction Kit** should be stored at room temperature in a cool, dry and well-ventilated area. All components of the kit can be adequately stored for up to 12 months.

The kit should be used in a well-ventilated area, keep away from the source of heat, sparks, open flames, and smoking.

To avoid evaporation, the pre-filled reagent should be used immediately after opening, and should not be placed for a long period of time.

Avoid exposure to UV light (e.g., for decontamination), which may result in accelerated aging.

## Sample Handling and Storage

Avoid foam inside or on the samples. Depending on the starting material, sample pre-treatment may be

required. Samples should be stored at room temperature (15~25°C) before starting the experiment. Samples should be used immediately after collection to extract nucleic acid or stored at 2~8°C for further experiment within 24 hours. For long-term storage, the samples should be placed at -20°C.

## Operation Guide

### 1. Reagent Preparation

**Proteinase K Solution Preparation:** Add 700µL Proteinase K diluent into Proteinase K, and make it completely dissolved, which can be used. Dissolved solution should be stored at -20°C or not more than 6 hours at room temperature, avoid repeated freezing and thawing (no more than 5 times).

### 2. Extraction Steps

- a. Add 20µL Proteinase K solution, 500µL Lysis solution and 200µL samples into 1.5mL nuclease-free centrifuge tube then vortex mixing for 15s. (Note: please do not use the solution which combines Proteinase K solution with Lysis solution).
- b. 60°C water bath for 20 mins, mixing several times during this time.
- c. Place the solution at room temperature for 3-5 mins, Instantly centrifuge to remove residual liquid from the top lid after the solution temperature in the same as room.
- d. Add 300µL absolute ethyl alcohol (self-contained) into centrifuge tube, and Instantly centrifuge to remove residual liquid from the top lid after mixing completely.
- e. Pipette 510µL solution from the above centrifuge, then add that into adsorption column; cover the lid and centrifuge at 8000rpm for 1 min; discard the waste liquid and put adsorption column into the collection tube.
- f. Repeat No.e once.
- g. Open the tube lid slowly, add 700µL washing A, centrifuge at 12000rpm for 1 min, then discard the waste liquid and put adsorption column into the collection tube.
- h. Open the tube lid slowly, add 700µL washing B, centrifuge at 12000rpm for 1 min, then discard the waste liquid and put adsorption column into the collection tube, after that centrifuge at 12000rpm for 2 mins.
- i. After centrifugation, the adsorption column was placed in a new 1.5mL nuclease-free centrifuge tube, then open the tube lid and place them at room temperature for 3-5 mins.
- j. Add 50µL eluent (should heating to 60°C) into adsorption, placed that at room temperature for 2-5 mins, centrifuge at 12000rpm for 1 min; then collected solution to 1.5mL centrifuge tube (self-contained), and stored it at -20°C.

## Troubleshooting Guide

This troubleshooting guide should assist you in resolving any problems that arise during the experimental process. For more information, please visit our Technical Support Centre and Frequently Asked Questions, page at: <http://www.medtl.net>. The scientists in our Tianlong company's Technical Services Department are always available to answer any questions you may have about the information and protocols contained in the manual, as well as sample and assay technologies (for contact information is included on the back cover or at: <http://www.medtl.net>).

When an exception or error occurs during the experiment, the current run step is terminated/stopped. After resolving the error or exception, restart the run from the beginning. The troubleshooting guide is shown in the following table.

| No. | Fault Symptom               | Fault Cause                         | Handling Method   |
|-----|-----------------------------|-------------------------------------|---|
| 1   | Poor extraction performance | Sample pretreatment                 | Please follow the operation requirements in the manual. |
|     |                             | The sample lysis may be incomplete. | Please follow the operation requirements in the manual. |
|     |                             | Other                               | Contact the after-sales service of our company.         |

\* Ensure that the reagents have been preserved and used according to the manufacturer's instructions.

## Quality Control

In accordance with Tianlong Company's ISO-certified Quality Management, each lot of **Viral DNA and RNA Extraction Kit** is tested against predetermined specifications to ensure consistent product quality.

## Limitations of Test Methods













The system performance has been established through performance evaluation studies using serum , swab etc samples to isolation viral DNA/RNA .

The user's responsibility is to validate system performance for any procedures performed in their laboratory that are not covered by Xi'an Tianlong Science and Technology Co., Ltd., performance evaluation studies.

Although the kit is intended for use in public health and scientific research, the purity and quality of extraction results are also affected by the testing instruments and personnel. Moreover, the kit uses a specially formulated eluent that can affect the absorbance value, so it is not recommended to use a UV rays spectrophotometer to measure the extraction effect directly.

The extraction kit is intended for use with clinical diagnostic samples, forensic materials, and scientific research samples. The instrument and operator have an effect on the concentration and purity of the extracted product. Any generated diagnostic results must be interpreted in conjunction with the other clinical or laboratory findings.

## Safety Symbols and Signs

| No. | Symbol  | Implication   |
|-----|---|---|
| 1   |    | Catalogue number                                    |
| 2   |    | Batch code  |
| 3   |    | Contains sufficient for <N> tests                   |
| 4   |    | Use by date   |
| 5   |    | Caution   |
| 6   |    | Temperature limit                                   |
| 7   |    | In vitro diagnostic medical device                  |
| 8   |    | Reminder  |
| 9   |    | Manufacturer  |
| 10  |    | Do not re-use                                       |
| 11  |  | Conformed with EU standard                          |
| 12  |  | Authorized representative in the European Community |

## Contact Information

For technical assistance and more information, please contact our Technical Support Center at +86-29-82210851 (Tel), +86-29-82216680 (Fax), [www.medtl.net](http://www.medtl.net) or contact your local distributor.

For up-to-date licensing information or product-specific disclaimers, please see the respective User Guide. Tianlong User Guides are available at [www.medtl.net](http://www.medtl.net) or can be requested from Tianlong Technical Services or the local distributor.

IFU\_T108H\_EN © 2021 Xi'an Tianlong Science and Technology Co., Ltd., all rights reserved.